T-SWAB TRANSPORT™

CTM-ADVANCE KIT

NUCLEIC ACID PRESERVATION & TRANSPORT MEDIUM FOR VIRUS (VIRUS INACTIVATION)

USER MANUAL FOR CTM-ADVANCE KIT

PRODUCT NAME CTM-Advance Kit

CAT. NO	제신 21-977호	제신 21-1009호
	PT2-NFS-1(1111) PT2-NFS-2(1112) PT2-NFS-3(1113) PT2-NFS-16(1115)	PT2-NFS-3G(1114)

1. Intended Use

CTM-Advance Kit consisting of the Swab and virus-inactivation medium is the specimen collection, nucleic acid preservation and transport medium kit for clinical virus such as influenza, SARS-COV-2 and etc. The preservation solution has the function of cell lysis and is compatible with the lysis buffer in the virus DNA/RNA extraction kit of all kinds of kits. In the subsequent nucleic acid extraction process, less or no lysate can be added.

2. Reagents

- Tris -Guanidine-HCL
- EDTA -N-Acetyl-L-Cysteine

3. Characteristics

The color of CTM-Advance Kit is colorless and each product contain 1 ml, 2 ml, or 3 ml of virus-inactivation media in the 12 ml vol. tube

4. NB Flocked Swab Type

NFS-1, NFS-2, NFS-3, NFS-16

5. Specimen Collection Procedures

- 1) After peel and open the pouch, collect specimen by using collection swab.
 - * Specimen collector has to hold the upper part of swab breakpoint, not lower part.
- 2) Collect appropriate amount of specimen coated on the surface of swab tip..
 - * The plastic applicator can be broken at breakpoint by excessive force.
- 3) After open the cap and put the swab specimen into the transport medium, carefully break the swab shaft against the breakpoint using both hands and discard the top portion of the swab shaft.
 - * Please be care careful when you break the breakpoint of swab without splash.
- 4) Recap the transport medium tube tightly and
- 5) Record the name of specimen, sex, age, number, date and etc.,
- 6) As soon as possible, transfer the tube containing specimen to the laboratory for analysis.
 - * If it is too late to deliver, store at 2-8°C and analyze in 72 hours. If stored at room temp., test it within 48 hours.

T-SWAB TRANSPORT™

CTM-ADVANCE KIT

NUCLEIC ACID PRESERVATION & TRANSPORT MEDIUM
FOR VIRUS (VIRUS INACTIVATION)

How to use CTM-Advance Kit

[Picture 1] Peel open the sealed pouch pack.

Remove one swab from the pouch and collect the specimen without bending the swab.

[Picture 2] Aseptically remove the cap from the vial and insert the swab into the tube containing medium.

Break the swab shaft by bending the swab.

[Picture 3] Replace the cap and secure the lid tight.

[Picture 4] Record the patient's information on the label.

Transfer the tube containing the specimen to laboratory for analysis.

6. Packing Unit

1 pc/individual package

50 pcs/Case

500 pcs/Carton Box

7. Storage & Expiration Date

Store: Room temp. (1 ~ 30°C)

Expiration Date: 24 months after production

8. Precautions

- Do not freeze or heat before use.
- Please don't use after expiry date.
- Do not use if it is damaged, broken, contaminated or leaking or the swab pouch is open.
- Prior to collecting please do not contact the reagent in tube.
- Single-use medical device.
- Do not reuse opened and used products.
- Record information on the label to deliver and analysis.
- Do not ingest the medium.
- Do not reuse, repackage, re-sterilize.

CREATED DATE: 2021.05.10

NOBLE BIOSCIENCES, INC.

13-50, SINBAEK-GIL, JEONGNAM-MYEON, HWASEONG-SI, GYEONGGI-DO, REPUBLIC OF KOREA
E-mail : info@noblebio.co.kr FAX : (82) - 31 - 291 - 0046 TEL : (82) - 31 - 291 - 0044 WWW.NOBLEBIO.CO.KR

NOBLE BIOSCIENCES, INC. IS THE COMPANY OF COLLECTION, TRANSPORT AND PRESERVATION SYSTEM FOR THE BIOPSY.

NASOPHARYNGEAL SWAB COLLECTION NASOPHARYNGEAL SPECIMEN COLLECTION METHOD


Information and Intended Use

Nasopharyngeal swabs are the preferred specimens for respiratory virus rapid antigen and RT-PCR testing. The nasopharyngeal swab has a white or semitransparent plastic shaft, with 3 different thicknesses, ending in a flock tip. There is a deep score mark (v) on the thick part as a breakpoint where it can be snapped into the transport medium. Each swab is individually packaged and labeled 'Noble Bio sterile flocked swab'.

Specimen Collection Procedures

1. Assemble all supplies such as powder-free gloves, mask, pen, facial tissues for patient, Nasopharyngeal flocked swab, transport medium and the other Lab requisition.
2. Check expiry date of transport medium and swab. And then, wash hands with soap and water or alcohol hand rub, and put on gloves and mask.
3. Have the patient sit in a chair or lie in on a bed and elevate the head of the bed so that their head can be tilted back up to 45°~75°. In the case of young children, please hold their head and chest, otherwise the swab shaft can be broken because of serious moving or sneezing.
4. Remove any mucous from the patient's nose with a tissue or cotton tipped swab prior to collecting specimens.
5. One half(1/2) of the distance from the corner of the nose to the front of the ear is usually about 4~5 cm(finest thickness of the swab shaft) to collect specimens.
6. Tilt the patient's head back slightly about 45°~75° to straighten the passage from the front of nose, and gently insert the swab along the base of the nose(not upwards) until it reaches the nasopharynx-gentle rotation of the swab may be helpful.(If resistance is encountered on one side, try the other nostril.) Rotate the swab several times to dislodge the columnar epithelial cells.
7. Put the swab specimens into the transport medium and break it at the breakpoint mark. The cap of the medium tube is screwed on tight.
8. Remove and discard gloves. Wash hands with soap and water or alcohol hand rub. Record all of the information about the specimens on the label. And send it to the laboratories according to their specifications.

T-SWAB TRANSPORT™

CTM-ADVANCE KIT

NUCLEIC ACID PRESERVATION & TRANSPORT MEDIUM FOR VIRUS (VIRUS INACTIVATION)

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Noble Bio
Clinical Diagnostic Products

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