
T-SWAB TRANSPORT™

CTM-ADVANCE SALIVA

NUCLEIC ACID PRESERVATION & TRANSPORT MEDIUM FOR VIRUS (VIRUS INACTIVATION)

How to use CTM-Advance SALIVA

[Picture 1] Peel open the sealed pouch pack.

Remove one swab from the pouch and collect the specimen without bending the swab.

[Picture 2] Refer to the **4. Specimen Collection Procedures** section of this manual, and collect saliva in the mouth.

[Picture 3] Aseptically remove the cap from the vial and insert the swab into the tube containing medium. Break the swab shaft by bending the swab.

[Picture 4] Replace the cap and secure the lid tight.

[Picture 5] Record the patient's information on the label.

Transfer the tube containing the specimen to laboratory for analysis.

5. Packing Unit

1 pc/individual package

50 pcs/Case

500 pcs/Carton Box

6. Reagents

- Tris
- Guanidine-HCL
- EDTA
- N-Acetyl-L-Cysteine
- Etc.

7. Storage & Expiration Date

Store: Room temp. (1 - 30°C)

Expiration Date: 24 months after production

8. Precautions

- Do not freeze or heat before use.
- Please don't use after expiry date.
- Do not use if it is damaged, broken, contaminated or leaking or the swab pouch is open.
- Prior to collecting please do not contact the reagent in tube.
- Single-use medical device.
- Do not reuse opened and used products.
- Record information on the label to deliver and analysis.
- Do not ingest the medium.
- Do not reuse, repack, re-sterilize.

CREATED DATE: 2021.08

T-SWAB TRANSPORT™

CTM-ADVANCE SALIVA

NUCLEIC ACID PRESERVATION & TRANSPORT MEDIUM FOR VIRUS (VIRUS INACTIVATION)

USER MANUAL FOR CTM-ADVANCE SALIVA

PRODUCT NAME CTM-Advance SALIVA

TYPES OF SWAB Saliva Swab

PERMIT NO. No. 21-1104

1. Intended Use

The CTM-Advance SALIVA consisting of Swab and virus inactivation medium is a sample collection, nucleic acid preservation and transport medium kit for clinical viruses such as influenza and SARS-CoV-2. The preservation solution is compatible with all kinds of viral DNA/RNA extraction kit lysis buffer with cell lysis function. There is no need to add additional lysis buffer in the subsequent nucleic acid extraction process. The types of Swab included in the kit is Saliva Swab with a cylindrical shape of plastic material for collecting saliva in the mouth.

2. Cat. No.

For Saliva Swab

PT2-SLS-1B(1122)

PT2-SLS-2B(1124)

PT2-SLS-1C(1123)

PT2-SLS-2C(1125)

Type of swab

Collection site

Saliva swab

Oral cavity

3. Characteristics

The color of CTM-Advance SALIVA medium is colorless and each product contain 2ml or 3 ml of virus-inactivation media in the 12 ml vol. tube

4. Specimen Collection Procedures

- 1) After peel and open the pouch, collect specimen by using collection swab.
 - * Specimen collector has to hold the upper part of swab breakpoint, not lower part.
- 2) Before collecting samples, massage both cheeks so that saliva is sufficiently filled in your mouth, and slowly swab through the stagnant saliva so that saliva can be sufficiently coated on the swab.
 - * The plastic applicator can be broken at breakpoint by excessive force.
- 3) After open the cap and put the swab specimen into the transport medium, carefully break the swab shaft against the breakpoint using both hands and discard the top portion of the swab shaft.
 - * Please be care careful when you break the breakpoint of swab without splash.
- 4) Recap the transport medium tube tightly and
- 5) Record the name of specimen, sex, age, number, date and etc.,
- 6) As soon as possible, transfer the tube containing specimen to the laboratory for analysis.
 - * If it is too late to deliver, store at 2-8°C and analyze in 72 hours. If stored at room temp., test it within 48 hours.


Noble Bio
Clinical Diagnostic Products

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Shape and Structure - Dimensions (Saliva swab used in the mouth)

Specimen collection tool (Permit No. 21-896)

No	Cat. No.	Full length A(mm)	Length from head to cut B(mm)	Head thickness C(mm)	Head length D(mm)	Stick thickness E(mm)	Handle thickness F(mm)
1	SLS-1 (1128)	150	43	7.5	28	2.6	2.6
2	SLS-2 (1129)	150	43	7.5	34	2.6	2.6

Margin of error: ±10%

Finished product picture (CTM-Advance SALIVA)



Component photos by model

No	Model No.	Pictures	Specimen transport medium/model
1	MAIN PT2-SLS-2B(1124)		2 ml / SLS-2(1129)
2	PT2-SLS-1B(1122)		2 ml / SLS-1(1128)
3	PT2-SLS-1C(1123)		3 ml / SLS-1(1128)
4	PT2-SLS-2C(1125)		3 ml / SLS-2(1129)

1. Massage both cheeks so that saliva is sufficiently filled in your mouth.

2. Slowly swab through the stagnant saliva so that saliva can be sufficiently coated on the swab.

NOBLE BIOSCIENCES, INC.

13-50, SINBAEK-GIL, JEONGNAM-MYEON, HWASEONG-SI, GYEONGGI-DO, REPUBLIC OF KOREA

E-mail : info@noblebio.co.kr

FAX : (82) - 31 - 291 - 0046

TEL : (82) - 31 - 291 - 0044

WWW.NOBLEBIO.CO.KR

NOBLE BIOSCIENCES, INC. IS THE COMPANY OF COLLECTION, TRANSPORT AND PRESERVATION SYSTEM FOR THE BIOPSY.

T-SWAB TRANSPORT™

CTM-ADVANCE SALIVA

NUCLEIC ACID PRESERVATION & TRANSPORT MEDIUM FOR VIRUS (VIRUS INACTIVATION)

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